

# 17 - 22PLUS

Simple solutions made to measure



# 17 - 22PLUS

## Technical data.

Technical data	17	22PLUS
Power supply voltage	220/240 V	220/240 V
Mains frequency	50/60 Hz	50/60 Hz
Rated power	2300 W (10 A)	2300 W (10 A)
External dimensions L x H x D	480 x 420 x 560 mm	480 x 420 x 660 mm
Total weight	55 kg	60 kg
Autonomy	From 6 to 10 cycles	From 5 to 9 cycles
Chamber dimensions Diam. x Depth	250 x 350 mm	250 x 450 mm

Cycle description	Cycle type	Holding time (min)	17	22PLUS
			Total cycle times (min) incl. drying	
134°C Universal	B	4	37	40
121°C Universal	B	20	52	55
134°C Prion	B	20	53	57
134°C Hollow	S	4	29	32
134°C Flash	S	4	25	27
XXX°C Custom	Temperature 121°C / 134°C Process times from 4' (134°C) or 20' (121°C), settable drying			
Vacuum test		TEST	27	25
Helix / B&D test		TEST	29	31

STANDARD EQUIPMENT	17	22PLUS
Usb port	•	•
Serial port (RS 232)	•	•
Bacteriological filter	•	•
Trays supplied	3	3
Tray extractor	•	•
Tray support	•	•

ACCESSORIES	17	22PLUS
Automatic filling kit	◦	◦
Auxiliary EV for <<demi>> water network	◦	◦
External printer	◦	◦

### Sede Legale ed Amministrativa / Head Quarter

CEFLA s.c.  
Via Selice Prov.le, 23/a  
40026 Imola (BO) - Italy  
Tel. +39 0542 653111  
Fax +39 0542 653344

### Fabbricatore / Manufacturer

Cefla Dental Group  
Via Bicocca, 14/C - 40026 Imola (BO) - Italy  
Tel. +39 0542 653441 - Fax +39 0542 653601  
www.sternweber.com  
sternweber@sternweber.com